

Plié 3 MPK Package - Order Form

Included in the package:

- Plié 3 knee, including 5 year warranty & the following items:
 - Wireless USB adapter
 - Lithium-ion battery case
 - Lithium-ion battery x2
 - Battery charger
 - Air pump with hose adapter
 - Smooth hose adapter for air pump
- Preferred designer fairing
- 1 x Freedom foot - Kinterra, Sierra, Highlander, Agilix, Pacifica LP or DynAdapt
- Set of adapters

How to order a Plié trial kit

Step 1	Step 2	Step 3	Step 4 - Successful Trial Only
Submit a completed order form to customerservices@steepergroup.com .	Once the PO and patient form has been received, this will generate an order for a 30 day trial period for the KIT-PLIEMPK. You will not be charged or invoiced at this stage.	Following a successful trial, call Customer Services on 0870 240 4133 to obtain an RTA number and return the demo Plié 3 knee only. If the trial was unsuccessful, call Customer Services on 0870 240 4133 to obtain an RTA number and return the full kit.	After contacting Customer Services your order will then be finalised. The definitive Plié 3 knee and the preferred designer fairing stated on this form will be delivered and invoiced.

N.B. If you do not require a trial and would like to place a definitive order, please submit this form, a purchase order using the part number KIT-PLIEMPK, and tick this box:

If only a foot is required:
If only a knee is required:

Patient Information:

Name:	Age:	Gender:	M	F	
Date knee is required for:					
Weight (kg) :	Height:	Skin tone:	Caucasian:	Mid Brown:	Brown:
Amputation Level (TF, TT, KD, HD...):	Left:	Right:	Size (22-31cm):	Bi-lateral:	
Impact level (see opposite for impact level descriptions):	Low:	Moderate:	High:		
Notes:					
Prosthetist:	Branch / Contractor:				

Plié 3 Knee Pyramid Adapter: Threaded Adapter:

Designer Fairings

Please select the preferred designer fairing. This will be supplied on the definitive Plié 3 knee upon completion of a successful evaluation and will not be shipped with the demonstration unit. The demonstration unit will be supplied in Standard Black.

- Standard Black
- Proud Marines
- Captivating Copper
- Radical Red
- The Patriot
- Pink Passion
- Forest Camo
- Fabulous Flames
- Tech Metal
- Blazing Blue
- Pearl White



Adapters

The following adapters are included within the MPK package. For alternative adapters, please place a separate order.

Set 1

Description	Part Number	Weight Rating	Material
Tube Clamp, 30 mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30 mm x 400 mm	NACC-00-12800-00	136kg	Titanium
Pyramid, 4 Hole, Female	NACC-00-13200-00	136kg	Titanium

Set 2

Description	Part Number	Weight Rating	Material
Tube Clamp, 30 mm, Female Pyramid	NACC-00-12500-00	136kg	Titanium
Tube w/Titanium Female Pyramid, 30mm x 400mm	NACC-00-12800-00	136kg	Titanium
Base, 3 Prong, Rotating Pyramid	NACC-00-13800-00	136kg	Stainless Steel
Pyramid, Rotating, Female	NACC-00-14000-00	136kg	Titanium

Freedom Foot Selection - select one



Category Selection: Kinterra			Category Selection - Sierra, Highlander, Agilix, Pacifica LP or DynAdapt			
Weight (kg)	Impact Level		Weight (kg)	Impact Level		
	LOW	MOD		LOW	MOD	HIGH
44-52	1	1	44-52	1	1	2
53-59	1	2	53-59	1	2	3
60-68	2	3	60-68	2	3	4
69-77	3	4	69-77	3	4	5
78-88	4	5	78-88	4	5	6
89-100	5	6	89-100	5	6	7
101-116	6	7	101-116	6	7	8
117-125	7		117-130	7	8	9
			131-147	8	9	
			148-166	9		

Side

Left: Right:

Toe Configuration¹

Regular Toe: Sandal Toe:

Size (cm)

22 23 24 25 26 27 28 29 30 31

Stiffness Category

1 2 3 4 5 6 7 8 9

Impact Level Descriptions

Low: Daily activities include mostly ground-level walking, moving around in the home and the community.

Moderate: Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

High: Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

¹Sandal toe available in sizes 22-28 for Kinterra, Sierra, Highlander, Agilix, Pacifica LP & DynAdapt Feet.

Plié 3 Warranty

5 years

The extended warranty on the Plié 3 from 3 to 5 years is automatically included within the MPK package. Please note, the Plié 3 knee must be returned for annual servicing to maintain the warranty. A loaner unit will be provided during this time.

For warranty terms on all other items, please refer to the Extended Warranty - KIT-PLIEMPK information.