


# evo+



Assistive Technology  
Environmental Control Device

 **steeper**  
Creating life's turning points, together



# evo+



evo+ is setting new standards for home automation with its ability to turn any iPhone, iPad or iPod Touch into a universal environmental control device.

Operated via touch screen, switches, or mice, evo+ is designed to provide independence to anyone who needs to remotely control their home devices.

The new design enables clients to personalise their home screen to access their most used icons, and adjust image sizes for hassle-free use. Seamless integration with other devices around the home, regular app updates and an improved intercom system, enables evo+ to provide a reliable and stable home control solution for any client.



## iOS compatible

Compatible with all iOS devices across iPhone, iPad and iPod ranges to offer continuous flexibility for clients.

## Seamless integration

evo+ integrates seamlessly with other peripheral devices around the home to create a full home control solution. Features such as improved intercom and dozens of app configuration options, make it a reliable partner in the client's home. See overleaf for the full compatibility of evo+.

## Reliable communication aid

evo+ can be used as a basic communication tool using frequently used, pre-recorded phrases.

## Simple & personalised setup

A Steeper engineer can deliver and install a basic evo+ system in just one hour, and will ensure that the app is configured according to the individual's specific daily requirements and preferences.

## Access options





evo+ can be accessed through touch screen or the use of a remote switch, so clients can navigate the bold icons easily and take control in their home.

## Customisable icons

Clients can select and prioritise icons based on their most used functions, ensuring they are clearly visible and easy to find. To further personalise the app, clients can pick their favourite colour theme from the full spectrum, and resize icons to increase their visibility and ease of access.

## Enhanced system

The enhanced evo+ app offers improved stability and reliability as it is fully integrated within the iOS system. Whenever a new iOS update is launched, clients will simply be prompted to update the evo+ app.



# What can I control with evo+?


evo+ provides clients full control of their devices, aided by our assistive technology engineers who create bespoke, automated solutions throughout their home. The ability to move around freely, communicate quickly, and adjust room conditions for comfort, allows individuals to spend more time embracing their independence and less time on daily tasks. See below the variety of peripheral devices that can be integrated with the evo+ app:



Wall/ceiling lights



Intercoms



Curtains/blinds



Profiling beds/  
chairs



Window openers



Automatic doors



Plug-in appliances  
e.g. lamps, fans, heaters



Home entertainment  
e.g. TV, DVD players



Telephones/pagers

A client's current assistive technology equipment can also be easily configured by our engineers to create a complete environmental control solution that runs throughout the home. The table below details the devices that can be used in conjunction with evo+, to create a safe and secure home individuals can maintain their independence within:

evo+		
	RF Bridge with Quad appliance	Minibirt
Wall/ceiling light		✓
Plug-in appliances e.g. fans, lamps, heaters		✓
Profiling bed/chair	✓	✓
Window opener	✓	✓
Curtain/blind controller	✓	✓
Automatic doors	✓	✓
Pager	✓	✓
Telephone		✓
Home entertainment e.g. TV, DVD player		✓
Intercom system		✓

For details on our full range of assistive technology devices or to request an assessment from the Steeper assistive technology team, please contact [assist@steepergroup.com](mailto:assist@steepergroup.com).