

## Arthritic Knee Sleeve



### Features and Benefits

- Simple knee sleeve
- Pull on design, elasticated for easy use of donning
- Provides compression and warmth for painful joints
- Posterior shaping eradicates bunching in popliteal area

### Indications for Use

- Arthritis

### Ordering Information

Size	Knee Circ. (cm)	Product Code
Small	31 - 34	803/AKS
Medium	34 - 37	805/AKS
Large	37 - 42	807/AKS
X-Large	42 - 46	808/AKS
XX-Large	49 - 54	809/AKS

## Stabilised Knee Sleeve



### Features and Benefits

- Incorporates medial and lateral flexible spiral stays for enhanced medial/lateral stability
- Provides compression and support
- Touch and close straps for improved fit
- Universal left/right application
- Made from Airprene

### Indications for Use

- Medio lateral stability and patella control
- knee instability
- Post-op rehabilitation
- Strains/sprains
- Bursitis
- Tendonitis osteoarthritis
- Rheumatoid arthritis
- Osteoarthritis

**Note:** When measuring the knee, ensure it is fully extended. For thigh circumference, measure 15cm above the centre of the patella. For calf circumference, measure 15cm below the centre of the patella.

### Ordering Information

Size	Thigh Circ. (cm)	Knee Circ. (cm)	Calf Circ. (cm)	Product Code
Small	34 - 36	25 - 28	26 - 28	SKS/CX/BE/S
Medium	36 - 39	28 - 30	28 - 32	SKS/CX/BE/M
Large	39 - 41	30 - 33	32 - 35	SKS/CX/BE/L
X-Large	41 - 44	33 - 36	35 - 37	SKS/CX/BE/XL

**More Information**  
Patient Information