

## Oasis TLSO




### Features and Benefits

Oasis is our most popular bivalve spinal system and ideal for patients who require maximum support, but cannot tolerate most rigid TLSO braces. Revolutionary spacer fabric offers increased breathability, heat and moisture removal, less bulk and longer life span than trilaminates. The Oasis TLSO includes a Dorsal Lumbar Kit (DLK) with Sternal Panel to prevent forward flexion. Effective between T7 - L5.

### Indications for Use

- Burst fracture
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Kyphosis
- Laminectomy
- Osteoporosis
- Paraplegia
- Postoperative
- Spinal stabilisation
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis



### Ordering Information

Size	Waist (cm)	Hip (cm)	Product Code
X-Small	66 - 82.3	78.7 - 94	OATL4-XS-PPK
Small	76.2 - 91.4	88.9 - 104.1	OATL4-S-PPK
Medium	86.4 - 101.6	99 - 114.3	OATL4-M-PPK
Large	96.5 - 111.8	109.2 - 124.5	OATL4-L-PPK
X-Large	106.7 - 121.9	119.4 - 134.6	OATL4-XL-PPK
XX-Large	116.8 - 132	129.5 - 144.8	OATL4-XXL-PPK
XXX-Large	127 - 142.2	139.7 - 154.9	OATL4-XXXL-PPK



#### More Information

Video: How to fit the Oasis TLSO brace

## Terrapin Bivalve TLSO Brace




### Features and Benefits

- The Bivalve technology allows for easy donning
- The Ilio-loc provides abdominal compression, positioning, and prevents migration
- Adjustable A/P control
- Removable and washable liner
- Made from a high density polyethylene
- Standard option has a front height of 10" (25.4cm)
- Short option has a front height of 8½" (21.6cm)
- Effective between T7 - L5

### Indications for Use

- Burst fracture
- Chronic back pain
- Compression fracture
- Degenerative, herniated or bulging disc
- Kyphosis
- Laminectomy
- Osteoporosis
- Paraplegia
- Postoperative
- Spinal stabilisation
- Spinal stenosis
- Spondylolysis
- Spondylolisthesis



### Ordering Information

Size	Hip Circ. (cm)	Hip Circ. (inch)	Product Code		
			Standard	Short	Pendulous Abdomen
1	76.2 - 86.4	30 - 34	TLSO-1	-	-
2	86.4 - 94	34 - 37	TLSO-2	-	-
3	94 - 104.1	37 - 41	TLSO-3	TLSO-3S	-
4	104.1 - 114.3	41 - 45	TLSO-4	TLSO-4S	TLSO-4P
5	114.3 - 124.5	45 - 49	TLSO-5	-	-



#### More Information

Video: How to fit the Terrapin TLSO brace