

Paediatric Hip Abduction Splint - Standard



- ### Features and Benefits
- Maintains concentric reduction of hips in patients with hip dysplasia
 - Replaces spica casts and harnesses for certain developmental dysplasia patients
 - One-Piece, non-adjustable polypropylene orthosis holds infant in 90° hip flexion and 60° abduction
 - Open design facilitates nappy changes
 - Closed-cell foam lining does not retain moisture

- ### Indications for Use
- Indicated for patients with hip dysplasia

Ordering Information

Size	Average Age (mnths)	Weight Guide (lbs)	Product Code
X-Small	<3	6 - 12	B/3899
Small	3 - 6	13 - 16	B/3900
Medium	6 - 12	17 - 21	B/3901
Large	12	22 - 30	B/3902
X-Large	18	30+	B/3903

Paediatric Hip Abduction Splint - Adjustable



- ### Features and Benefits
- Lightweight, paediatric orthosis offers advantage of abduction adjustment over the standard paediatric hip splint model described opposite
 - Sturdy polypropylene frame set in 90° hip flexion - abduction easily adjustable 40° through 180° by aluminium slotted screw assembly
 - Fits infants from 3 months (6-12lbs) to 18 months (over 30lbs)
 - More total contact than Pavlik harness
 - Easy application, removal facilitates nappy changes
 - Closed-cell foam lining does not retain moisture

- ### Indications for Use
- Indicated for patients with hip dysplasia
 - X-Large size effective for certain cerebral palsy patients

Ordering Information

Size	Average Age (mnths)	Weight Guide (lbs)	Product Code
X-Small	<3	6 - 12	B/3906
Small	3 - 6	13 - 16	B/3907
Medium	6 - 12	17 - 21	B/3908
Large	12	22 - 30	B/3909
X-Large	18	30+	B/3910