

Push Med Ankle Brace



Features and Benefits

- Provides long-term support of the ankle joint
- Based on the functions of tape bandage
- The degree of pressure and pronation/supination limits is adjustable
- Silicone application ensuring good position retention
- Comfort zones located on the most vulnerable areas of the ankle such as the Achilles tendon
- Made of Sympress™ which is a high quality material made for ultimate comfort
- Microfibres keep the skin dry and ensure optimal comfortability

Indications for Use

- Moderate residual instability
- Treatment of acute lateral ankle ligament lesions (alternative to taping)
- Treatment of stable ankle fractures
- Aftercare for ankle fractures
- Postoperative treatment of ligament injuries and fractures
- Capsular irritation and tendonitis following ankle fracture or ligament lesions
- Arthrosis with instability
- Secondary prevention of ankle ligament lesions (under moderate stress with accumulation of fluid)
- Impingement



Ordering Information

Size	Ankle Circ. (cm)	Product Code	
		Left	Right
1	26 - 29	PUS022B	PUS0017B
2	29 - 32	PUS023B	PUS0018B
3	32 - 35	PUS024B	PUS0019B
4	35 - 38	PUS025B	PUS020B
5	38 - 41	PUS026B	PUS021B

Push Aequi Ankle Brace



Features and Benefits

- An effective ankle orthosis which strikes a balance between efficient mechanical support of a disabled ankle and good user comfort
- The rigid support is situated exclusively on the medial side of the ankle and has a recess for the malleolus
- The lateral side is made from preformed foam, which easily adapts to individual anatomy

Indications for Use

- Serious (residual) instability (under high stress)
- Treatment of acute lateral ankle ligament lesions
- After-care of conservative or operated ankle fractures
- Arthrosis with instability without accumulation of fluid
- Secondary prevention of ankle ligament lesions

Note: Paediatric sizes also available, see **Page 228**.

Ordering Information

Size	Short Heel Circ. (cm)	Product Code	
		Left	Right
1	27 - 31	320111	320121
2	31 - 34	320112	320122
3	34 - 40	320113	320123