

ToeOFF® Short



Features and Benefits

- Specifically developed to accommodate individuals with shorter stature, slightly broader-shaped tibia, and wider foot at the metatarsal area
- ToeOFF® Short has a heel height that is 7mm in all sizes

Indications for Use

- Designed to manage footdrop in conditions such as
 - Stroke
 - Multiple Sclerosis
 - Post Polio Syndrome
 - Muscular Dystrophy
 - Spinal Cord Injuries
 - Traumatic Brain Injuries
 - Guillian-Barre Syndrome
 - Charcot-Marie-Tooth
 - Myelomeningocele
 - Neuropathy
 - Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations



Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Footplate Length (cm)	Height (cm)	Product Code	
					Left	Right
Small	1 - 2.5	32 - 35	21.5	29.5	289901011	289902011
Medium	2.5 - 5	35 - 38	23.5	32	289901012	289902012
Large	5 - 8	38 - 42	25.5	34	289901013	289902013

ToeOFF® 2.0



Features and Benefits

- The ToeOFF® 2.0 is the updated version of the ToeOFF® Original
- Lightweight, carbon reinforced rehabilitation orthosis that absorbs force on heel strike
- Worn at the anterior side of the leg, bridging the ankle and preventing plantar flexion
- SoftKIT™ 2.0 comprises padding shaped to precisely fit the orthosis and comes supplied with two tibia pads
- ComfortKIT™ 2.0 comprises thicker padding, designed for the more sensitive skin, fits the orthosis precisely
- Durable and strong straps
- Available in black

Indications for Use

- Mild to severe footdrop
- Mild to moderate ankle instability
- Impaired knee control
- Limb proprioception deficit
- Mild to moderate spasticity
- Partial foot amputation



*This accessory comes in a pair.

Note: To order ToeOFF® Original please contact Customer Services.

Ordering Information

Size	Shoe Size (UK)	Footplate Length (cm)	Height (cm)	Product Code	
				Left	Right
X-Small	1 - 3	21	36	289201010	289202010
Small	3 - 4	23	38	289201011	289202011
Medium	5 - 9	24.5	40.5	289201012	289202012
Large	10 - 11	27	43	289201013	289202013
X-Large	12 - 13	28.5	45	289201014	289202014

Accessories			
Size	Product Code		
	SoftKIT™ 2.0*	ComfortKIT™ 2.0*	Replacement Straps
X-Small	287500010	287510010	287600010
Small	287500011	287510011	287600011
Medium	287500012	287510012	287600012
Large	287500013	287510013	287600013
X-Large	287500014	287510014	287600014

ToeOFF® 2½



Features and Benefits

- Features the same dynamic functions and stability as ToeOFF® 2.0 but has a lower heel height of 7mm, to adapt to modern shoes
- Lower toe lift which results in more space in the toe box
- MikroFIX pre-applied to simplify the padding application
- Durable and strong straps
- The straps can be applied on the interior or external side of the AFO

Indications for Use

- Designed to manage footdrop in conditions such as:
 - Stroke
 - Multiple Sclerosis
 - Post Polio Syndrome
 - Muscular Dystrophy
 - Spinal Cord Injuries
 - Traumatic Brain Injuries
 - Guillian-Barre Syndrome
 - Charcot-Marie-Tooth
 - Myelomeningocele
 - Neuropathy
 - Cerebral Palsy
- Supports gait in conditions such as Posterior Tibialis Tendon Dysfunction (PTTD) and toe amputations
- Can also be used for partial foot amputations, most proximal level is Chopart

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Footplate Length (cm)	Height (cm)	Product Code	
					Left	Right
X-Small	1 - 2.5	32 - 35	21.5	36	289801010	289802010
Small	2.5 - 5	35 - 38	23	38	289801011	289802011
Medium	5 - 8	38 - 42	24.5	40.5	289801012	289802012
Large	8 - 10.5	42 - 45	27	43	289801013	289802013
X-Large	> 0.5	>45	28.5	45	289801014	289802014

ToeOFF® Flow 2½



Features and Benefits

- ½ heel height of 7mm in all sizes
- Smoother transition (flow throughout the gait cycle)
- Increased ROM in sagittal plane
- Accommodates lower shoe heel heights
- Accommodates more gentle contours of shoe insoles
- More clearance for forefoot in the shoe toe box
- Easier to customise to meet alignment and biomechanical needs
- MikroFIX pre-applied which makes it quicker to attach a SoftKIT™ 2.0, ComfortKIT™ 2.0 or CoverKIT™ 2.0
- Closure straps are the same as for ToeOFF® 2.0
- Lightweight and thin
- Open heel



Indications for Use

- Mild proximal deficits
- Mild spasticity
- Limb proprioception deficit
- Functions well for higher activity levels

Ordering Information

Size	Shoe Size (UK)	Shoe Size (EU)	Footplate Length (cm)	Height (cm)	Product Code	
					Left	Right
X-Small	1 - 2.5	32 - 35	21.5	36	289601010	289602010
Small	2.5 - 5	35 - 38	23	38	289601011	289602011
Medium	5 - 8	38 - 42	24.5	40.5	289601012	289602012
Large	8 - 10.5	42 - 45	27	43	289601013	289602013
X-Large	>10.5	>45	28.5	45	289601014	289602014