

	Name	Impact Level	K Level	Features and Benefits	Foot Size (cm)	Build Height on size 25cm	Weight based on size 25cm inc foot shell	Weight Limit	Warranty	Part Number	Trial Period	Sandal Toe Available
Everyday	Agilix™		K3	<ul style="list-style-type: none"> <li>Vertical shock absorption and carbon fibre flexibility optimised for low to moderate impact K3 ambulators; reduces socket shear forces and improves overall comfort</li> <li>Split keel design for increased durability</li> <li>The low profile and lightweight design of the Agilix allows users to do more and feel less tired at the end of the day</li> </ul>	22-31	14.3cm	553g	166kg	3 Years	NF15 ● ▲ ■	30 days (only sizes 23-30cm, categories 2-8)	Yes*
	Sierra®		K3	<ul style="list-style-type: none"> <li>Full length heel ensures exceptionally smooth gait</li> <li>Split keel design</li> <li>Integral male pyramid</li> <li>Available in 9 stiffness categories</li> </ul>	22-31	14.6cm	515g	166kg	3 Years	NFS1 ● ▲ ■	30 days (only sizes 23-30cm, categories 2-8)	Yes*
	Highlander®		K3	<ul style="list-style-type: none"> <li>Low build height with integral male pyramid</li> <li>Split keel design for excellent inversion/eversion</li> <li>Weight limit 166kg</li> <li>Available in 9 stiffness categories</li> </ul>	22-31	10.8cm	470g	166kg	3 Years	NFS3 ● ▲ ■	30 days (only sizes 23-30cm, categories 2-8)	Yes*
	Highlander® Max		K3	<ul style="list-style-type: none"> <li>Low build height with integral male pyramid</li> <li>Split keel design for excellent inversion/eversion</li> <li>Weight limit 226kg</li> <li>Available in 3 stiffness categories</li> </ul>	26-31	12.1cm (26cm)	590g (26cm)	226kg	1 Year	NFS3 ● ▲ ■ (categories 10-12)	No trial period	Yes*
	DynAdapt™		K3-4	<ul style="list-style-type: none"> <li>Slim profile, carbon fibre foot system with a slender, anatomical design</li> <li>Multi-axial function for maximum comfort</li> <li>EnduraCore Technology delivers up to 3 times the fatigue life of a standard laminate</li> <li>Full length split keel gives an uninterrupted rollover with good compliance and a dynamic response through to toe-off</li> </ul>	22-31	15.8cm	535g	166kg	3 Years	NF10 ● ▲ ■	30 days (only sizes 23-30cm, categories 2-8)	Yes*
	Shockwave		K3-4	<ul style="list-style-type: none"> <li>Spherical Shock Unit (SSU) provides ultimate vertical shock absorption at heel strike and also a transverse rotation of +/- 15 degrees easing stress and pressure on the residual limb</li> <li>Full-length heel initiates plantar flexion at heel strike for better ground compliance and maximised energy return</li> </ul>	23-30	17.6cm	897g	147kg	3 Years	NF20 ● ▲ ■	Please speak to a Steeper Product Manager to organise a trial	Yes*
Low-Profile Carbon Foot	Pacifica LP®		K3	<ul style="list-style-type: none"> <li>Ultra low build height</li> <li>Integral male pyramid</li> <li>Split keel design</li> <li>Foot lies within the foot shell to improve cosmetic finishing</li> <li>Available in 9 stiffness categories</li> </ul>	22-31	4.7cm	465g	166kg	3 Years	NFS4 ● ▲ ■	30 day trial (only sizes 23-30cm, categories 2-8)	Yes*
Hydraulic Ankles	Kinterra		K2-3	<ul style="list-style-type: none"> <li>Articulating ankle delivers reduced socket pressures for healthier residual limbs</li> <li>Increased ground contact and a more symmetrical gait</li> <li>Improved stability and safety, especially on slopes and varied terrain</li> <li>Dorsi-Assist spring aids ground clearance in swing phase</li> <li>Individual adjustment of dorsi and plantar flexion</li> </ul>	22-31	12cm	795g	125kg	3 Years	NROM ● ▲ ■	30 day trial (only sizes 23-30cm, categories 2-7)	Yes*
Microprocessor Ankle/Foot	Kinnex 2.0		K2-3	<ul style="list-style-type: none"> <li>Sensor-controlled ankle movement reads and processes data at a rate of 200 times per second and Rapid Response™ value technology responds at a rate of 5 milliseconds</li> <li>IP67 rated - up to 30 minutes in 1 metre of fresh water</li> <li>New hydraulic oil reservoir and deal technology for increased oil capacity and longevity</li> </ul>	24-30	18.4cm	1432g	125kg	3-5 Years	NF14 ● ▲ ■	Please speak to a Steeper Product Manager to organise a trial	Yes*

## Spare Parts Ordering Information

Spectra® Sock	Colour	Size (cm)	Sock Size (cm)
NSO-NPS-20030-00	B	22-25	30
NSO-NPS-20035-00	B	26-28	35
NSO-NPS-20040-00	B	29-31	40

Colour: B = Black

Note: These part numbers apply to all PROTEOR Feet except for Kinnex 2.0.

Kinnex 2.0		
Spectra® Sock	Size (cm)	Sock Size (cm)
NSO-F14-00024-00	24	24
NSO-F14-00025-00	25	25
NSO-F14-00026-00	26	26
NSO-F14-00027-00	27	27
NSO-F14-00028-00	28	28
NSO-F14-00029-00	29	29
NSO-F14-00030-00	30	30

Pacifica LP	
Bumper Kit	Size (cm)
NKIT-00-11323-00	22-25
NKIT-00-11324-00	26-28
NKIT-00-11325-00	29-31

Highlander/Highlander Max	
Bumper Kit	Size (cm)
NKIT-00-11303-00*	22-25
NKIT-00-11304-00	26-28
NKIT-00-11305-00	29-31

\*This size is not available for the Highlander Max.

Agilix	
Bumper Kit	Size (cm)
NKIT-00-11363-00	22-25
NKIT-00-11364-00	26-28
NKIT-00-11365-00	29-31

Sierra	
Bumper Kit	Size (cm)
NKIT-00-11353-00	22-25
NKIT-00-11354-00	26-28
NKIT-00-11355-00	29-31

DynAdapt	
Bumper Kit	Size (cm)
NKIT-00-11373-00	22-25
NKIT-00-11374-00	26-28
NKIT-00-11375-00	29-31

Shockwave	
Bumper Kit	Size (cm)
NKIT-00-11463-00	23-25
NKIT-00-11464-00	26-28
NKIT-00-11465-00	29-30

Kinterra	
Bumper Kit	Size (cm)
NKIT-00-11333-00	22-25
NKIT-00-11334-00	26-28
NKIT-00-11335-00	29-31

Kinnex 2.0	
Bumper Kit	Size (cm)
NKIT-00-11403-00	24-25
NKIT-00-11404-00	26-28
NKIT-00-11405-00	29-30

Foot Shell	
Part Number	Foot Shell Type
NFTC-3M-0 ▲ ● 4-R ■	Standard, with CAP
NFTC-2M-1 ▲ ● 4-S ■	Sandal Toe, without CAP

Note: These part numbers apply to all PROTEOR feet except for Shockwave and Kinnex 2.0.

Foot Shell - Shockwave	
Part Number	Foot Shell Type
NFTC-3F-0 ▲ ● 4-R ■	Standard, with CAP
NFTC-2F-1 ▲ ● 4-S ■	Sandal Toe, without CAP

Foot Shell - Kinnex 2.0	
Part Number	Foot Shell Type
NFTC-3K-1 ▲ ● 4-R ■	Standard, with CAP
NFTC-2K-1 ▲ ● 4-S ■	Sandal Toe, without CAP

Spare Parts Ordering Guide  
 ▲ = Size in cm (22-31 = Regular, 22-28 = Sandal Toe)  
 ● = Skin Tone (L = Light, M = Medium, D = Dark)  
 ■ = "L" or "R" for Left or Right

Spare Parts Ordering Guide - Shockwave  
 ▲ = Size in cm (23-30 = Regular, 23-28 = Sandal Toe)  
 ● = Skin Tone (L = Light, M = Medium, D = Dark)  
 ■ = "L" or "R" for Left or Right

Spare Parts Ordering Guide - Kinnex 2.0  
 ▲ = Size in cm (24-30 = Regular, 24-28 = Sandal Toe)  
 ● = Skin Tone (L = Light, M = Medium, D = Dark)  
 ■ = "L" or "R" for Left or Right

## Impact Level Descriptions

**Low:** Daily activities include mostly level ground walking, moving around in the home and the community.

**Moderate:** Daily activities include up to unlimited walking, climbing stairs and occasional moderate lifting (this does not include running).

**High:** Daily activities include fast walking, jogging, running, lifting heavy objects and/or recreational sports.

Note: Although the 'high' impact level states 'running' as an appropriate activity, the Pacifica LP is not ideal for running.

= **Water resistant**  
Safe for contact with water and occasional submersion. Rinse immediately with fresh water and dry following exposure to water, particularly salt or chlorinated.

= **Waterproof**  
The Kinnex should not be immersed in salt water or chlorinated water - if exposed to salt or chlorinated water, immediately rinse with fresh water and dry.

Note: For sandal toe option add an "S" before the side of foot e.g. when ordering Kinterra category 1, size 22 sandal toe left foot, the part number would be NROM122SL

**Order Guide**

● = Category (1-12)\*\*  
 ▲ = Size in cm (22-31)  
 ■ = "L" or "R" for Left or Right

\*\*Category range differs for each foot, please refer to the Steeper Lower Limb Catalogue for category tables.

\*Sandal toe option only available in sizes 22-28. Please refer to the Steeper Lower Limb Catalogue for individual sandal toe options.