

## Elegance Plus Gloves

We have adapted and developed our Elegance silicone gloves, and are pleased to introduce the new Elegance Plus cosmetics to our cosmetic glove range. These enhanced gloves are designed to fit the Steeper Myo Kinisi myoelectric hand, Steeper cable operated hands, and functional hands from other manufacturers.

Crafted using multiple layers of advanced silicone material and constructed with an integrated, reinforced inner mesh layer, these gloves are highly tear-resistant. This silicone material allows for additional stretch in comparison to a PVC glove, and so reduces the energy expended by the myoelectric hand motor each time it opens - making this a highly robust cosmetics particularly suited to higher activity patients.

For training and more details on how to order visit [www.upskillbysteepers.com](http://www.upskillbysteepers.com)




### Features and Benefits

- New integrated, reinforced inner mesh layer provides increased resilience
- Enhanced cosmetic finish with additional detailing on the palms, knuckles, nails and joints
- Improved tear resistance and highly robust
- New easy-glide coating for ease of use when putting on clothes
- Compatible with Steeper and other manufacturers' myoelectric and mechanical hands
- Exclusively available in the advanced micropigmentation TrueFinish™ silicone for added realism and depth
- Increased longevity, making this range highly cost effective

### Elegance Plus Colour Options

\*Only available in TrueFinish

E0	E1	E2	E3	E4	E5	E6	E7	E8	E9	E10	E11	E12	E13	E14	E15	E16	E17	E18
----	----	----	----	----	----	----	----	----	----	-----	-----	-----	-----	-----	-----	-----	-----	-----



## Elegance Plus Ordering Information

If a non-standard finish or bilateral glove is required, please state at the time of ordering.

**Note:** # Denotes colour.

Glove Part Number	Left/Right	Glove Number	Palm Circumference	Length
TSG306/E#-EP	Left	306	184mm (7 ¼)	645mm (25 ¼)
TSG307/E#-EP	Right	307	184mm (7 ¼)	645mm (25 ¼)
TSG304/E#-EP	Left	304	184mm (7 ¼)	645mm (25 ¼)
TSG305/E#-EP	Right	305	184mm (7 ¼)	645mm (25 ¼)
TSG312/E#-EP	Left	312	197mm (7 ¼)	620mm (24 ¼)
TSG313/E#-EP	Right	313	197mm (7 ¼)	620mm (24 ¼)
TSG302/E#-EP	Left	302	197mm (7 ¼)	645mm (25 ¼)
TSG303/E#-EP	Right	303	197mm (7 ¼)	645mm (25 ¼)
TSG300/E#-EP	Left	300	210mm (8 ¼)	645mm (25 ¼)
TSG301/E#-EP	Right	301	210mm (8 ¼)	645mm (25 ¼)

Swatches are available free of charge

### Ordering Examples

TSG300/E13-EP - Orders an 8 ¼ Left Silicone Elegance Plus Glove in TrueFinish™ colour E13

TSG305/E7-EP - Orders and 7 ¼ Right Silicone Elegance Plus Glove in TrueFinish™ colour E7