



medipro®
New product Information

medipro® catalogue supplement Stump Socks medipro® RELAX Night Care

The optimum supplement to the medipro® Liner RELAX and RELAX TF/ TFC



Allivates phantom pain, sensations and idiopathic residual limb pain without side effects!

INDICATIONS

- Protection of the residual limb
- Phantom pain and phantom sensations
- Idiopathic residual limb pain

NEW: Now available in two different lengths!

PROPERTIES

- Relieves or eliminates phantom pain, phantom sensations and idiopathic residual limb pain due to isolation from electromagnetic effects thanks to the metallurgical knit (Umbrellan®)
- Double-knit: protects the residual limb against cold, light bumps, friction, other external influences
- Knit with temperature-regulating thread (Clima-Comfort®):
 - warming and cooling function, depending on environmental temperature
 - highly elastic knit with SOFT topbands for comfortable, reliable fit – without slipping or constriction
 - pleasant wearing comfort due to the soft, skin-friendly Umbrellan® and Clima-Comfort® material
 - no heat build-up, as stocking is breathable and moisture-regulating thanks to the well-known Clima-Comfort® technology
 - particularly suitable for wearing at night and any time no prosthesis is being worn

TECHNICAL DATA

Material	Umbrellan® with Clima-Comfort®
Sizes / Length	see product table
Amputation level	Below-knee, knee disarticulation/above-knee
Cleaning	Hand wash at max. 30°C using medipro® care neutral soap M. 870.902
Supplied	Umbrellan® User Set (medipro® pain diary, medipro® pain scale)



Stick onto Page 140 of the medipro® General Catalogue



Umbrellan® is a unique, patented knitted fabric that combines maximum elasticity with an excellent electromagnetic screen.

The materials used on the metallurgical base are hypoallergenic, climate-regulating and have antibacterial properties.

Its functionality, therapeutic benefit and quality have been statistically significantly demonstrated.

Prospective randomised study available as Art. No. 70.914.



clima-comfort®



140a

140b

TEXTILES VERTRAUEN
Schadstoffgeprüfte Textilien
Tested for harmful substances
nach / according to
Öko-Tex Standard 100



medipro®
New product information


medipro® catalogue supplement

Stump Socks

medipro® RELAX Night Care



The optimum supplement to the medipro® Liner RELAX and RELAX TF/ TFC



NEW! Now even quicker and easier to measure for size!

SIZE SELECTION

(Below-knee/knee exarticulation/above-knee)

1. Choose the suitable stocking length depending on the **length of the residual limb L** (30 or 40 cm)*
2. Measure the **circumference A** (at top band)
3. The residual limb length and circumference give the correct medipro® RELAX Night Care size

* Important note:

- medipro® RELAX Night Care is available in two standard lengths (30 or 40 cm). Please choose the nearest size in relation to the length of the residual limb.
- The stocking can be carefully sewn up for residual limbs < 30 cm.



Art. No. M.870.902

medipro® care neutral soap, 200 ml

For thorough cleansing of the residual limb and the stocking.

Best if ordered together with the stocking.

99.802/09.07/3,5'



4 102639814308 12

PRODUCT TABLE

Designation	Circumference A (level of topband)	Size	Art. No. Length 30 cm	Art. No. Length 40 cm	Unit	Material
RELAX Night Care	20 – 26 cm	S	–	1.880.322	Piece	Umbrellan®
RELAX Night Care	26 – 33 cm	M	1.880.328 NEW!	1.880.323	Piece	Umbrellan®
RELAX Night Care	33 – 42 cm	L	1.880.329 NEW!	1.880.324	Piece	Umbrellan®
RELAX Night Care	42 – 52 cm	XL	1.880.325	1.880.327 NEW!	Piece	Umbrellan®
RELAX Night Care	52 – 62 cm	XXL	1.880.326	–	Piece	Umbrellan®

Maximum useful life 3 months.