

KEO

Take control of your independence

An easy to use, wireless, environmental remote control system designed to facilitate independent living for people who lack full mobility.

Lightweight and portable, the KEO can be configured according to an individual's unique needs and enables a huge range of everyday tasks such as opening curtains, answering telephones, unlocking doors and changing TV channels, giving a new kind of freedom to the user in the home environment.

Flexibility:

The KEO can be programmed to hold up to approximately 250 infrared codes to control many different appliances including:

- TV, video, stereo and other appliances with an infrared remote control
- Loud speaking telephone using infrared control
- Intercoms
- Mobile phones – Using Keo with Bluetooth
- Computers – Using the IR Receiver for PC



KEO



Options available:

KEO - infrared output only.

KEO with Bluetooth - infrared and Bluetooth output. Allows mobile phone control.

IR Receiver for PC - Use with either KEO, KEO with Bluetooth or ACTIV 500. The IR Receiver plugs into the USB port of your computer and allows you to use your KEO to control your computer.

KEO Switch Transmitter - Plug a switch into the Switch Transmitter and access KEO or KEO with Bluetooth using switch scanning.

Features:

- Infrared and Bluetooth control
- Light and portable
- Touch screen and switch access

Assistive Technology Division,
Unit 6, Hunslet Trading Estate,
Severn Road, Leeds LS10 1BL

Tel: +44 (0) 113 207 0449,
Fax: +44 (0) 113 270 9346
Email: assist@rslsteeper.com
www.rslsteeper.com