



Local details



Passive hand and cosmetic glove range

Prosthetic product solutions



Passive hand and cosmetic glove range

We are pleased to introduce the next generation of passive hands and matching cosmetic gloves. Our advanced materials and manufacturing make these detailed, durable and lightweight passive hands and cosmetic gloves.



Passive hand features

- Now 20 percent lighter.
- Robust construction for easy digit positioning.
- Three built-in wrist options for easy connection to oval or round wrists, or 12mm threaded studs.
- Five adult shapes and sizes.
- Available with extended wires and without internal moulding, to make prostheses for partial hands and wrist disarticulation.

Cosmetic glove features

- All-new moulding detail.
- 560mm – 620mm sleeve length.
- Silicone or PVC.
- Available in 19 colours and in TrueFinish™ or Standard ranges.

The hands are suitable for the following connectors:

Round wrist	
Hand plate to fit wrist B01155	B01172
Male stud 1/2" x 20 TPI	B26431
Male stud M12 x 1.5mm	B26430
Friction wrist hand plate	B12817
Fixed wrist connector: laminated / glued	B05439

Oval wrist	
Male stud 1/2" x 20 TPI	B26431
Male stud M12 x 1.5mm	B26430
Endoskeletal, fixed	B06856
Fixed wrist connector: laminated / glued	B06255

A bracket for socket or palm case lamination when using the passive hand with extended wires is available:

Part number	B05864
-------------	--------

Ordering information for passive hands

Size UK	Length (mm)	Wrist dia. (mm)	Application	Left part no.	Right part no.	Ext. wire left part no.	Ext. wire right part no.
3 1/2	140	54	Male	FHL1	FHL2	FHEXL1	FHEXL2
3 1/4	140	54	Male	FHL3	FHL4	FHEXL3	FHEXL4
3 1/4	140	54	Male long fingers	FHL21	FHL22	FHEXL21	FHEXL22
3	140	54	Male/Female	FHL5	FHL6	FHEXL5	FHEXL6
2 3/4	140	54	Female	FHL7	FHL8	FHEXL7	FHEXL8

Extended wire and standard passive hands are both suitable for the same cosmetic glove as denoted by their part number.

Silicone TrueFinish™ gloves for passive hands

UK size	Approx. length (mm)	Application	Glove left part no.	Glove right part no.	Colour range part no.	Hand left part no.	Hand right part no.	Euro size (mm) palm circumference	Euro size (mm) middle finger length
3 1/2	600	Male	TSGL1	TSGL2	E0 to E18	FHL1	FHL2	212	100
3 1/4	620	Male	TSGL3	TSGL4	E0 to E18	FHL3	FHL4	205	89
3 1/4	620	Male long fingers	TSGL21	TSGL22	E0 to E18	FHL21	FHL22	202	95
3	600	Male/Female	TSGL5	TSGL6	E0 to E18	FHL5	FHL6	186	83
2 3/4	560	Female	TSGL7	TSGL8	E0 to E18	FHL7	FHL8	172	82

To order standard silicone glove colours, omit the T prefix. These are available for the range E0 to E18. For example: SGL2/E7

Ordering example: passive hand FHL2 + TSGL2/E2. This provides a 3 1/2 right passive hand with TrueFinish™ 3 1/2 right silicone glove in colour E2.

PVC TrueFinish™ gloves for passive hands

UK size	Approx. length (mm)	Application	Glove left part no.	Glove right part no.	Colour range part no.	Hand left part no.	Hand right part no.	Euro size (mm) palm circumference	Euro size (mm) middle finger length
3 1/2	600	Male	TGL1	TGL2	E0 to E18	FHL1	FHL2	212	100
3 1/4	620	Male	TGL3	TGL4	E0 to E18	FHL3	FHL4	205	89
3 1/4	620	Male long fingers	TGL21	TGL22	E0 to E18	FHL21	FHL22	202	95
3	600	Male/Female	TGL5	TGL6	E0 to E18	FHL5	FHL6	186	83
2 3/4	560	Female	TGL7	TGL8	E0 to E18	FHL7	FHL8	172	82

To order standard PVC glove colours, change part number prefix to CGL. These are available for the range E0 to E18. For example: CGL8/E4

Ordering example: passive hand FHL1 + TGL/E5. This provides a 3 1/2 left passive hand with a TrueFinish™ 3 1/2 left PVC glove in colour E5.

Switch ordering information E24654: E0 to E18 Standard colours
E24655: E0 to E18 TrueFinish™ colours