

REPLACEMENT OF EXPLORER WITH LOW ACTIVITY LINER

Due to the supplier of the outer fabric of Explorer products (pin and cushion liners, and suspension sleeves) going out of business, Silipos has been unable to source the exact same material.

For this reason we have taken the opportunity to source a new material, which will have the side benefit of increasing the durability of the new LA liner.

The Low Activity Liner is the replacement product for the Explorer. It provides the same skin conditioning benefits and comforts as Explorer, as it uses the exact same gel. The LA liner is also manufactured on the exact same moulds; therefore there should be no issues when replacing an Explorer for an LA liner.

The table below shows the part numbers for Explorer Pin liners, and the corresponding LA liner part number that should be ordered going forward.

Explorer Part number	Size/Style	LA liner part number
N18001	Medium 3mm Uniform	N18101
N18002	Medium 3mm Tapered	N18102
N18003	Medium 6mm Uniform	N18103
N18004	Medium 6mm Tapered	N18104
N18005	Medium 9mm Uniform	N18105
N18006	Medium 9mm Tapered	N18106
N18007	Medium + 3mm Uniform	N18107
N18008	Medium + 3mm Tapered	N18108
N18009	Medium + 6mm Uniform	N18109
N18010	Medium + 6mm Tapered	N18110
N18011	Medium + 9mm Uniform	N18111
N18012	Medium + 9mm Tapered	N18112
N18013	Large 3mm Uniform	N18113
N18014	Large 6mm Uniform	N18114
N18015	Large 9mm Uniform	N18115
N18016	Large + 3mm Uniform	N18116
N18017	Large + 6mm Uniform	N18117
N18018	Large + 9mm Uniform	N18118
N18019	X large 3mm Uniform	N18119
N18020	X Large 6mm Uniform	N18120


The table below shows the part numbers for Explorer Cushion liners, and the corresponding LA liner part number that should be ordered going forward.

Explorer Part Number	Size/Style	LA liner part number
N18041	Medium 3mm Uniform	N18141
N18042	Medium 3mm Tapered	N18142
N18043	Medium 6mm Uniform	N18143
N18044	Medium 6mm Tapered	N18144
N18045	Medium 9mm Uniform	N18145
N18046	Medium 9mm Tapered	N18146
N18047	Medium + 3mm Uniform	N18147
N18048	Medium + 3mm Tapered	N18148
N18049	Medium + 6mm Uniform	N18149
N18050	Medium + 6mm Tapered	N18150
N18051	Medium + 9mm Uniform	N18151
N18052	Medium + 9mm Tapered	N18152
N18053	Large 3mm Uniform	N18153
N18054	Large 6mm Uniform	N18154
N18055	Large 9mm Uniform	N18155
N18056	Large + 3mm Uniform	N18156
N18057	Large + 6mm Uniform	N18157
N18058	Large + 9mm Uniform	N18158
N18059	X Large 3mm Uniform	N18159
N18060	X Large 6mm Uniform	N18160



I attach the sizing chart for measuring for a new LA liner, but again please understand that if you are currently using an Explorer you can swap it out for an LA liner.

Sizing Guidelines



Liner Size	Proximal Circumference (Measure 15" from Distal End)			
	Proximal		Distal	
	Relaxed	Fully Stretched	Relaxed	Fully Stretched
Medium	9" / 23cm	16" / 41cm	8" / 20cm	12" / 30cm
Medium +	13" / 33cm	21" / 53cm	9" / 23cm	12" / 30cm
Large	15" / 38cm	22" / 56cm	9" / 23cm	15" / 38cm
Large +	16" / 41cm	27" / 68cm	12" / 30cm	16" / 41cm
X Large	16" / 41cm	30" / 76cm	12" / 30cm	20" / 51cm

Distal Circumference (Measure 1.5" from Distal End)

The Explorer Suspension sleeves will also be completely discontinued, and will be superseded by the DuraGel Suspension sleeves. The table below gives the current Explorer sleeve item numbers, with the corresponding DuraGel sleeve number which should be ordered going forward.

There were only 3 sizes in Explorer sleeves, but we are introducing a 4th size for the DuraGel range, Large.

Regular (un-reinforced) Sleeves

Explorer	Size	DuraGel
N18081	Medium	N18181
N18082	Medium +	N18182
No item number	Large	N18183
N18083	Large +	N18184

Reinforced Sleeves

Explorer	Size	DuraGel
N18091	Medium	N18191
N18092	Medium +	N18192
No item number	Large	N18193
N18093	Large +	N18194



Tel. +44 (0) 870 240 4133
 Fax. +44 (0) 870 240 4779
 Email. sales@rslsteeper.co.uk
 Web. www.rslsteeper.com